

A more efficient, disposable articulator

A new option to send completed bites to your clients.

by Vic Sober

Durable and stronger than most other plastic articulators, Wondertech Articulators (from Dental Creations) do not have stability problems. The articulators function much like metal "Foster" articulators. Once articulated, the case can easily be put through the same

functions and excursive movements as metal or semi-adjustable articulators.

Many other plastic articulators have a tendency to rock back and forth or they can flex and bow, which makes it difficult to work. The dental technician constantly has to

reposition the articulator to make it stable. Wondertech Articulators are made of a more durable, stable plastic, which helps to make a better end product.

The Wondertech Articulator functions the same as all other popular small, fixed metal articulators



Wondertech Small Articulators

Wondertech Articulators

A disposable and affordable alternative to metal articulators.

Features include:

- Affordable & disposable
- Light weight & easily adjustable
- Made of durable ABS plastic
- Provide stable mounting

Dental Creations

wonderfill.com | 254-772-4661

CIRCLE RS #18




Fig. 1 The components to the Wondertech Articulators



Fig. 2 Apply small amount of adhesive to model base



Fig. 3 Center upper ball receiver attachment




Fig. 4 Assemble ball joint and articular base together. Snap ball joint into ball receiver.



Fig. 5 Upper model assembled




Fig. 6 Slide lower arm assembly into articulator base




Fig. 7 Place bite or triple tray onto model




Fig. 8 Set bite




Fig. 9 Place lower assembly into place. The articulator will move around through the ball joint to aid in articulation.




Fig. 10 Glue lower assembly onto lower model



Fig. 11 Glue upper ball joint together



Fig. 12 Glue lower assembly onto articulator base



Fig. 13 Wondertech Articulator fully assembled

and is disposable, which eliminates the messy clean-up that is required with the metal articulators. Of remarkable importance, the Wondertech Articulators are affordable enough that the technician can send the finished product to the dental office with the articulator attached. They are also lighter weight than a metal articulator, which saves in delivery costs. By sending an articulator with the finished product, if there is a problem, the dentist can see where the mistake occurred. Metal articulators are too expensive, so they cannot be sent to the dentist with the finished product attached.

Additionally, Wondertech Articulators are designed to be easier to work with because whether a lab technician is doing crown and bridge work or a full denture, the Wondertech Articulators can be easily used for both applications. When working with metal articulators, the technician either has to use the metal crown and bridge articulator or use the full denture articulator to make the application.

Finally, Wondertech Articulators are easier to mount. The bottom piece slides in and out and the top piece will move up and down or from left to right which makes it easier to adjust the model to the right angle. With these articulators, the technician can send the finished product to the dentist with the articulator attached. After simply trimming the excess plaster from the outside of the product, articulator and product can be sent to the dentist.

With metal articulators—and some plastic options—before sending the finished work to the dentist, it must be removed from the articulator, the excess plaster knocked off and the model trimmed while ensuring the model is not broken or cracked. The Wondertech Articulators save bench time and, therefore, money for the dental lab.

Below is a step-by-step technique for using Wondertech Articulators, a great option for any dental technician looking to find an affordable and efficient alternative to traditional, metal articulators.

STEP-BY-STEP: USING THE ARTICULATOR

01 Using a cyanoacrylate gap filling medium glue, apply the adhesive to the model base. This will hold the model in place for proper viewing of the final case. (**Fig. 2**)

02 To prevent a misaligned articulation, center the upper ball receiver attachment to the middle of the model base. This will ensure a successful result when it's time to send the case to the client. (**Fig. 3**)

03 “Snapping” in the ball joint into the ball receiver attachment assures you have properly connected the two. Your upper model is correctly assembled because the models will be level with the bench. (**Figs 4 & 5**)

04 Assemble the lower arm assembly by aligning and sliding the two arms with the two holes in the articulator base. (**Fig. 6**)

05 Trim bite or triple tray and place onto the models. This step makes certain the final model will match the bite. (**Fig. 7**)

06 Set the bite by bringing the two models together with the bite in place. (**Fig. 8**)

07 Now place the lower assembly into place in its place below the prepared upper portion. Via the flexible ball joint, the articulator will move around to aid in articulation. (**Fig. 9**)

08 Using a cyanoacrylate gap filling medium glue, adhere the lower assembly onto the lower model. (**Fig. 10**)

09 Once you know the bite is set, you may glue the upper ball joint assembly into place. (**Fig. 11**)

10 To finish the articulation, glue the two arms of the lower assembly to the two holes of the articulator base. (**Fig. 12**)

Using this method, lab technicians can ensure excellent results with Wondertech Articulators. This affordable, disposable solution can help any lab ensure they deliver reliable, esthetic results to their clients every time. **lab**

ABOUT THE AUTHOR

Vic Sober is a technical manager at a lab in Texas and has been a crown and bridge technician for nearly 20 years. He attended the Naval Dental School of Dental Assisting and Technology in San Diego, Cali. and graduated in 1985. Throughout the years Mr. Sober has been a PFM team leader, Dental Laboratory Operations Manager and a retired Navy Chief Petty Officer. He continues to bring a high degree of excellence to the dental industry.