

WONDER Z-RING KIT

Wonder Z-Ring creates an even heat distribution across the firing tray which provides for perfectly glazed restorations regardless where they are placed on the tray. Kit includes: 1 adjustable radiant heat ring, 1 honeycomb firing tray with 20 stainless steel pins, and 3 fiber free firing trays.

ADVANTAGES:

- Eliminates cold zones on firing tray
- Completely even heat across the tray surface
- Eliminates dull glaze finishes (fire scale) caused by uneven temperature in firing zone
- Optimizes depth of shade vitality of all layered ceramic applications
- Adjustable ring accommodates larger cases

Item # 101 Retail Price \$69.95



1.254.772.4661 / www.wonderfill.com
[YouTube.com/DentalCreationsLtd](https://www.youtube.com/DentalCreationsLtd)