

WONDERTECH SILICONE BASE FORMERS

Wondertech Silicone Base Formers are made from silicone rubber. Easy to use – simply pour stone, level off the stone, and then place the teeth cast on the form. They are reusable. They come in a set of 5 per box – Small Full Arch, Medium Full Arch, Large Full Arch, Right and Left Quadrants.

Item No. 430
Retail Price \$27.95

 *from the makers of Wonderfill®*
dentalcreations ltd.
1.254.772.4661 / www.wonderfill.com
[YouTube.com/DentalCreationsLtd](https://www.youtube.com/DentalCreationsLtd)