

SASSY SILICONE MIXING CUPS

Sassy Silicone Mixing Cups are acrylic non-stick cups made of a special silicone rubber surface adaptable for dental acrylics, resins and cements. The cups are ideal for premixing acrylic because their wide base prevents spillage and allows for instant removal of congealed acrylic. They are resistant to high temperatures.

The mixing cups are now available in three sizes: Small (1.25" Dia. x 1" H, holds 5 cc), Medium (2" Dia. x 1 ½" H, holds 30 cc), Large (2 ¾" Dia. x 2 ⅕" H, holds 80 cc).

Sassy Silicone Mixing Cups Product Line

Small (Set of 2)	DC Item No. 774	Suggested Retail Price	\$8.50
Medium	DC Item No. 775	Suggested Retail Price	\$7.50
Large	DC Item No. 776	Suggested Retail Price	\$11.95

