


# Seeing the light on optimized shade matching

How the Shade Wand helps eliminate ambient light variables to help match shades and minimize remakes. [compiled by Noah Levine]

*Information provided by Dental Creations.*



## THE PROBLEM

**Environmental factors make matching a patient's tooth shade a challenge**

Depending on the light in the room, the clinician's vision and even the viewing angle the clinician is using, tooth shades can appear different, making it difficult to make a precise match.

## THE SOLUTION

### Shade Wand

- The 5,500 kelvin LED light generates neutral, ambient light at three power settings
- Provides the optimal amount of light at the patient's mouth for easier and more precise shade matching
- Eliminates color distraction and outside color sources from interfering with shade selection
- Powered by rechargeable batteries that provide up to 4 hours of full power light time on a single charge
- No warm-up time required before it's ready to use

Dental Creations Ltd.  
254-772-4661 | wonderfill.com  
CIRCLE RS #XX