

SMART DAPPEN DISH
Premium Glass Mixing Dish

Smart Dappen Dish is the easy and effective way to mix acrylics, amalgam, resins, sealants, and other dense materials. This premium glass mixing dish is specially designed with a smooth, transparent, heat resistant glass to prevent materials from sticking to the dish. Featuring a thick solid base, this dappen dish offers high stability to protect from spills. All glass dappen dishes are non-autoclavable, cold sterilization only. These dishes are offered in a variety of colors: amber, clear, blue, and green, but each dish is sold individually.

ADVANTAGES:

- Mix acrylics, amalgam, resins, sealants, and other dense materials
- Smooth, transparent, heat resistant glass prevents sticking
- Thick base prevents spilling
- Cold sterilization cleaning only
- Available colors: amber, clear, blue, & green
- Sold individually

Clear Smart Dappen Dish

Item #765

Suggested Retail Price \$3.50

Amber Smart Dappen Dish

Item #766

Suggested Retail Price \$3.75

Blue Smart Dappen Dish

Item #767

Suggested Retail Price \$3.75

Green Smart Dappen Dish

Item #768

Suggested Retail Price \$3.75

