

## ***WONDERWHEELZ***

### ***High Quality Silicone Rubber Wheels!***

Wonderwheelz square edge silicone rubber wheels are specially formulated with the finest grade of silicone and are excellent for pre-polishing porcelains and metals. Wonderwheelz are very flexible and adapt to contours. Because Wonderwheelz are long lasting, they are cost effective.

#### **Item# 501 (Blue Silicone Wheel Box of 20)**

Fine – for high luster polishing Size: 22 x 3.15 mm / 7/8" x 1/8"

Suggested Retail Price \$16.49

#### **Item# 511 (Blue Silicone Wheel Box of 100)**

Fine – for high luster polishing Size: 22 x 3.15 mm / 7/8" x 1/8"

Suggested Retail Price \$50.95

#### **Item# 502 (Black Silicone Wheel Box of 20)**

Medium – for smoothing Size: 22 x 3.15 mm / 7/8" x 1/8"

Suggested Retail Price \$16.49

#### **Item# 512 (Black Silicone Wheel Box of 100)**

Medium – for smoothing Size: 22 x 3.15 mm / 7/8" x 1/8"

Suggested Retail Price \$50.95

#### **Item# 503 (White Silicone Wheel Item Box of 20)**

Course – for reducing and smoothing Size: 22 x 3.15 mm / 7/8" x 1/8"

Suggested Retail Price \$16.49

#### **Item# 513 (White Silicone Wheel Item Box of 100)**

Course – for reducing and smoothing Size: 22 x 3.15 mm / 7/8" x 1/8"

Suggested Retail Price \$50.95

