

***PERFECT PEGS***  
***Ceramic Firing Pegs***

Perfect Pegs, ceramic firing pegs, hold crowns and bridges in place in a porcelain furnace. They provide stability and allow for close placement. Perfect Pegs are excellent for use with Emax and Zirconia restorations.

***ADVANTAGES:***

- Excellent for use with Emax and Zirconia restorations
- Allow close placement for crowns & bridges
- Provide stability
- Can be fired up to 2200° F
- Come in 3 different shapes

**Item #112** (For Small Anteriors / Small Bicuspid Pack of 6)  
Suggested Retail Price \$26.95

**Item #114** (For Small Molars Pack of 6)  
Suggested Retail Price \$26.95

**Item #116** (For Large Molars / Large Bicuspid Pack of 6)  
Suggested Retail Price \$26.95

**Item #118** (For Posterior Pack of 6)  
Suggested Retail Price \$26.95

**Item #119** (Assortment | Pack of 12)

- 4 Small Anterior / Small Bicuspid pegs
- 4 Small Molars pegs
- 4 Large Molar / Large Bicuspid pegs

Suggested Retail Price \$45.55

