

PERFECT PEGS
Ceramic Firing Pegs

Perfect Pegs, ceramic firing pegs, hold crowns and bridges in place in a porcelain furnace. They provide stability and allow for close placement. Perfect Pegs are excellent for use with Emax and Zirconia restorations.

ADVANTAGES:

- Excellent for use with Emax and Zirconia restorations
- Allow close placement for crowns & bridges
- Provide stability
- Can be fired up to 2200° F
- Come in 3 different shapes

Item #112 (For Small Anteriors / Small Bicuspids Pack of 6)
Suggested Retail Price \$26.95

Item #114 (For Small Molars Pack of 6)
Suggested Retail Price \$26.95

Item #116 (For Large Molars / Large Bicuspids Pack of 6)
Suggested Retail Price \$26.95

Item #118 (For Posteriors Pack of 6)
Suggested Retail Price \$26.95

Item #119 (Assortment | Pack of 12)

- 4 Small Anterior / Small Bicuspids pegs
- 4 Small Molars pegs
- 4 Large Molar / Large Bicuspids pegs

Suggested Retail Price \$45.55

