

***SASSY FLEXIBLE MIXING BOWLS***  
***For A Bubble Free Mix***

Sassy Flexible Mixing Bowls allow for an easy, more efficient, and bubble-free mixing of alginate, plaster, stone, and investments. These bowls are resistant to most chemicals, solvents and oils making them easy to clean! Just flex the bowl and then wipe it clean! This allergy approved, latex free flexible bowl is a must have addition to any office or lab. Sassy Flexible Mixing Bowls can be used with most automatic alginate mixers.

***Advantages:***

- Flexible design allows for inside-out folding
- Latex free
- Bowls can be used with most automatic alginate mixers
- More flexible than rubber bowls
- Self-cleaning, just flex & wipe off
- Autoclavable option available

Autoclavable Black – Medium – 450 CC  
Item #780  
Suggested Retail Price \$7.75

Non-Autoclavable Purple – Medium – 500 CC  
Item #782  
Suggested Retail Price \$7.15

Non-Autoclavable Black – Medium – 500 CC  
Item #783  
Suggested Retail Price \$7.15

