

WONDERFILL

Tongue & Void Filler

The pre-mixed EZ in, EZ out Tongue and Void Filler eliminates time spent mixing alginate filler or grinding out stone from lower impressions. It also fills in bubbles, air pockets, missing teeth, and undercuts to enhance the fit of appliances and bleaching trays. It can be used as a temporary dam when making partial models. It holds up under heat when making heat-shrink retainers.

ADVANTAGES:

- Eliminates time spent mixing alginate filler or grinding out stone from lower impressions
- Fills in bubbles, missing teeth, & undercuts to enhance the fit of appliances & bleaching trays
- Use as a temporary dam when making partial models
- Holds up under heat when making heat-shrink retainers

Item# 344

Suggested Retail Price \$54.95

