

## ***WONDERWHEELZ***

### ***Diamond Impregnated Rubber Wheels***

Wonderwheelz Diamond Impregnated Rubber Wheels are excellent for cutting, abrading, removing scratches, and smoothing zirconia, porcelain, lithium disilicate, metals, and hard materials.

Recommended speed is 8,000-12,000 rpm.

#### **Item# 515 (Red Diamond Impregnated Rubber Wheel)**

Coarse Grade – for cutting and abrading zirconia, porcelain, lithium disilicate, metals, and hard materials.

Wheel dimensions: 26mm x 2mm 7/8"

Suggested Retail Price \$49.95

#### **Item# 516 (Blue Diamond Impregnated Rubber Wheel)**

Medium Grade – for removing deep scratches and smoothing the surface of zirconia, porcelain, lithium disilicate, metals, and hard materials.

Wheel dimensions: 26mm x 2mm 7/8"

Suggested Retail Price \$49.95

