

WONDERWHEELZ

High Quality Silicone Rubber Wheels!

Wonderwheelz square edge silicone rubber wheels are specially formulated with the finest grade of silicone and are excellent for pre-polishing porcelains and metals. Wonderwheelz are very flexible and adapt to contours. Because Wonderwheelz are long lasting, they are cost effective.

Item# 501 (Blue Silicone Wheel Box of 20)

Fine – for high luster polishing Size: 22 x 3.15 mm / 7/8" x 1/8"

Suggested Retail Price \$17.95

Item# 511 (Blue Silicone Wheel Box of 100)

Fine – for high luster polishing Size: 22 x 3.15 mm / 7/8" x 1/8"

Suggested Retail Price \$52.95

Item# 502 (Black Silicone Wheel Box of 20)

Medium – for smoothing Size: 22 x 3.15 mm / 7/8" x 1/8"

Suggested Retail Price \$17.95

Item# 512 (Black Silicone Wheel Box of 100)

Medium – for smoothing Size: 22 x 3.15 mm / 7/8" x 1/8"

Suggested Retail Price \$52.95

Item# 503 (White Silicone Wheel Item Box of 20)

Course – for reducing and smoothing Size: 22 x 3.15 mm / 7/8" x 1/8"

Suggested Retail Price \$17.95

Item# 513 (White Silicone Wheel Item Box of 100)

Course – for reducing and smoothing Size: 22 x 3.15 mm / 7/8" x 1/8"

Suggested Retail Price \$52.95

