

PERFECT CHEEK RETRACTORS

For A Clear View

Perfect Cheek Retractors are the new handsfree, self-retaining way to get a full view of your patient's mouth during direct bonding of anterior and posterior teeth. This break resistant, clear plastic allows for excellent intraoral and extraoral pictures by providing an unobstructed view of the oral cavity. These Cheek Retractors are compatible with cold sterilization or autoclavable up to 275° F. Sold 2 per pack.

ADVANTAGES:

- Break resistant, clear plastic
- Handsfree & self-retaining
- Autoclavable up to 275° F or cold sterilization
- Helpful for intraoral or extraoral pictures
- Available in adult & child size

Perfect Cheek Retractors – Autoclavable – Pack of 2 - Adult Item #810

Suggested Retail \$16.95

Perfect Cheek Retractors - Autoclavable – Pack of 2 - Child Item #812

Suggested Retail \$16.95

