

***SHADE WAND***  
***LED Shade Matching Light***

The Shade Wand is a hand held, rechargeable LED shade matching light that brings the perfect light to the patient's mouth. This 9-inch Shade Wand uses a 5500 kelvin LED light source which generates the proper amount of light at the patient's mouth. This enables the doctor to perfectly match tooth shades. The Shade Wand offers three different light settings and is constructed with quality components to ensure long life. It gives four full hours of light time at 100 percent light and eight full hours at 50 percent light.

(9" X 1")

***ADVANTAGES:***

- Rechargeable, cordless & easy to use
- 5500 kelvin LED light with 3 light settings
- Eliminates color distraction in oral cavity caused by outside color sources interfering with shade selection
- No warm up time required
- Pays for itself in one remake

Item #163  
Suggested Retail Price \$129.99

