

WONDERPINZ

#2 Medium Brass Stick Dowel Pins

- Non-corrosive, high quality brass stick pins
- Ideal for all impression material
- Model pin for easy fixation in the impression
- Long spike ensures secure hold in the impression
- Easily slide to end position via the conical pin shape
- Easy to remove the spike without damaging the die
- Great for production of anatomical and saw cut models
- Plaster arch and model base can be made one after the other in the impression

Item #730 (#2 Medium Brass Stick Pins – 500 pieces)
Suggested Retail Price \$59.95

Item #735 (#2 Medium Brass Stick Pins – 100 pieces)
Suggested Retail Price \$16.95

