

Z-BEADS

Zirconia Sintering Beads

Z-Beads Zirconia Sintering Beads provide a better temperature rate when using sintering trays. Because the sintering tray and Z-Beads heat at two different temperature rates, Z-Beads absorb the heat and provide a controlled temperature which helps prevent shrinkage. Zirconia Sintering Beads support crowns, inlays, onlays, veneers, and bridges during the sintering process. Z-Beads can withstand up to 1700 degrees Celsius and are made of high purity zirconia that provide contamination free results for restorations. For best results, use with Vented Sintering Tray.

Advantages:

- Helps control temperature of sintering tray
- Perfectly sized at 2mm
- Contamination free firing results
- Controls shrinkage
- High purity
- Withstands up to 1700 degrees Celsius

Size: 200g Jar of 2mm Beads

Item #146

Suggested Retail Price: \$45.50

