

## ***ALUMINUM OXIDE***

Aluminum Oxide is a dental grade abrasive that comes in 45 micron, 280 grit, that can be used for micro abrasion in dental applications and used for divesting and cleaning appliances. It is a fine grade, non-contaminating pure aluminum oxide for use in pressure blasting systems and air brushes. It is 99.42% (white) pure electro-fused Aluminum Oxide. This product is only available for purchase in the U.S.A.

### ***ADVANTAGES:***

- Easy to use
- Saves time
- Rough surface preparation
- Smoothing/polishing
- Bulk removal of investment and tough oxide layers
- Alloy preparation for porcelain work

#### **Aluminum Oxide - 1lb.**

**Item#: 850**

**Suggested Retail Price: \$37.95**

#### **Aluminum Oxide - 2 lbs.**

**Item#: 851**

**Suggested Retail Price: \$59.95**

