

## ***PERFECT CHEEK RETRACTORS***

### ***For A Clear View***

Perfect Cheek Retractors are the new handsfree, self-retaining way to get a full view of your patient's mouth during direct bonding of anterior and posterior teeth. This break resistant, clear plastic allows for excellent intraoral and extraoral pictures by providing an unobstructed view of the oral cavity. These Cheek Retractors are compatible with cold sterilization or autoclavable up to 275° F. Sold 2 per pack.

### ***ADVANTAGES:***

- Break resistant, clear plastic
- Handsfree & self-retaining
- Autoclavable up to 275° F or cold sterilization
- Helpful for intraoral or extraoral pictures
- Available in adult & child size

### **Perfect Cheek Retractors – Autoclavable – Pack of 2 - Adult Item #810**

Suggested Retail \$17.95

### **Perfect Cheek Retractors - Autoclavable – Pack of 2 - Child Item #812**

Suggested Retail \$17.95

