

***PERFECT PEGS – Small Anteriors / Small Bicuspids
Ceramic Firing Pegs***

Perfect Pegs, ceramic firing pegs, hold crowns and bridges in place in a porcelain furnace. They provide stability and allow for close placement. Perfect Pegs are excellent for use with Emax and Zirconia restorations.

ADVANTAGES:

- Excellent for use with Emax and Zirconia restorations
- Allow close placement for crowns & bridges
- Provide stability
- Can be fired up to 2200° F
- Come in 3 different shapes

Item #112 (For Small Anteriors / Small Bicuspids Pack of 6)
Suggested Retail Price \$27.50

