

***SMART DAPPEN DISH***  
***Premium Glass Mixing Dish***

Smart Dappen Dish is the easy and effective way to mix acrylics, amalgam, resins, sealants, and other dense materials. This premium glass mixing dish is specially designed with a smooth, transparent, heat resistant glass to prevent materials from sticking to the dish. Featuring a thick solid base, this dappen dish offers high stability to protect from spills. All glass dappen dishes are non-autoclavable, cold sterilization only. These dishes are offered in a variety of colors: amber, clear, blue, and green, but each dish is sold individually.

***ADVANTAGES:***

- Mix acrylics, amalgam, resins, sealants, and other dense materials
- Smooth, transparent, heat resistant glass prevents sticking
- Thick base prevents spilling
- Cold sterilization cleaning only
- Available colors: amber, clear, blue, & green
- Sold individually

**Clear Smart Dappen Dish**

Item #765

Suggested Retail Price \$4.25

**Amber Smart Dappen Dish**

Item #766

Suggested Retail Price \$4.25

**Blue Smart Dappen Dish**

Item #767

Suggested Retail Price \$4.25

**Green Smart Dappen Dish**

Item #768

Suggested Retail Price \$4.25

