

***PERFECT PEGS***  
***Ceramic Firing Pegs***

Perfect Pegs, ceramic firing pegs, hold crowns and bridges in place in a porcelain furnace. They provide stability and allow for close placement. Perfect Pegs are excellent for use with Emax and Zirconia restorations.

***ADVANTAGES:***

- Excellent for use with Emax and Zirconia restorations
- Allow close placement for crowns & bridges
- Provide stability
- Can be fired up to 2200° F
- Come in 3 different shapes

**Item #119 (Assortment | Pack of 12)**

- 4 Small Anterior / Small Bicuspids pegs
- 4 Small Molars pegs
- 4 Large Molar / Large Bicuspids pegs

Suggested Retail Price \$41.95

