

ZIR-HOLDER

Crown Holder with Silicone Tips

Zir-Holder is used to hold crowns easily and more conveniently. Zir-Holder has soft silicone tips that will protect crowns from scratches. This crown holder is light and comfortable to use. The Zir-Holder handle has multi-purpose functionality. It comes with a ball carver, a scoop, and a millimeter ruler. The handle also comes with a vibration tip for brushes.

ADVANTAGES:

- 2-in-1 Anterior & Posterior crown holder
- Silicone tips protects from scratches
- Handle includes scoop, ball carver and ruler
- Easily & conveniently holds crowns in place

Handle Size: 4.15 x 0.35 in (105 x 8.90 mm)

Item# 161

Suggested Retail Price \$65.95

