

SMART SALIVA EJECTORS

These disposable saliva ejectors are built for dependable performance and patient comfort. The flexible tubing includes an internal wire that allows the ejector to bend easily and stay in the desired position without collapsing. A soft, non-removable tip is securely bonded to the tube, ensuring it remains firmly in place during use. The side-opening design provides effective suction while helping protect oral tissues from being drawn in. Available in clear with a blue tip, each bag contains 100 ejectors.

ADVANTAGES:

- Secure, bonded tip
- Soft and comfortable
- Consistent suction
- Holds shape when bent

Smart Saliva Ejectors – Pack of 100

Item #: 42

Suggested Retail Price: \$7.95

